


## Creating cocktails for <sup>19</sup>F NMR-assisted FBDD approach based on Enamine's in-stock collection

D. Volochnyuk, K. Melnykov, Ya. Filatov, S. Tyukhtenko, S. Ryabukhin

## Status quo of <sup>19</sup>F NMR-assisted FBDD



## TAXS (FLUORINE CHEMICAL SHIFT ANISOTROPY AND EXCHANGE FOR SCREENING)

- One of the biophysically available and highly effective tools for analyzing the presence of ligand-protein
- Fast and requires only a limited amount of protein, and therefore compares favorably with other established non-NMR methods used in high throughput screening
- Although <sup>1</sup>H detection was used in the original competition-based approaches, there are significant advantages of using <sup>19</sup>F detection
- The absence of spectrum overlap allows you to check large chemical mixtures and automatically analyze spectra even in the presence of protonated buffers, solvents, and detergents.
- The problems with human/semiautomated data mining (especially for diastereomeric "CF<sub>2</sub>" cpd)
- Doubling of signals in amide spectra (hindered rotation) and NH-pyrazoles (slow tautomeric exchange)
- Doubling of signals due to their instability in the aquatic environment (oxadiazoles)
- Solubility problems for some compounds
- The presence of a mixture of diastereomers and compounds with two non-equivalent fluorines
- Not everyone passed the filters for unwanted compounds

## Cocktails for <sup>19</sup>F NMR Fragment Screening: From Theory to Application

