




Safe, Multigram Flow Generation of Diazomethane for Organic Synthesis

Serhiy Ryabukhin, Viacheslav Pendiukh, Alexander Ye. Pashenko, Dmytro Volochnyuk


Background of the project

Dry Diazomethane in Flow: Safety, Large Scale, Wide Diversity of Products

- A flow system for generating diazomethane enables the synthesis of MedChem building blocks and intermediates with productivity up to 0.45 mol/h, scalable to hundreds of grams per operation.
 - **Diazoketones**, synthesized using this system, are valuable, bench-stable reagents¹ on their own, and can be used for **synthesis of 2,4-thiazoles² and pyrazoles**, as well as other cyclization reactions.
 - The use of flow photoreactor enables the synthesis of **β-amino acids³**, **azetidinones**, **piperidinones**, and **cyclopropenes**, with the **preservation of chirality** of the starting acids for diazoketones and their derivatives.




Amino Acid-Derived Diazoketones (AADDKs) Thermal Stability Summary: Safe, Bench Stable




- AADDKs remain stable after a year at room temperature, showing promise as **bench-stable synthetic equivalents for halogenoketones**.
 - TGA of diazoketones show their decomposition between 110-142 °C, exothermically, indicating **medium thermal risk and reaction hazard indices**.

$T(1^{\text{st}} \text{ process}), ^\circ\text{C}$	137	142
$T(2^{\text{nd}} \text{ process}), ^\circ\text{C}$	160	167
$E_A(1^{\text{st}} \text{ process}), \text{kJ/mol}$	220	115
$T, ^\circ\text{C}$ range for EA	140-153	142-160

Flow-Generated Diazomethane in Reaction Cascade: Setup and Scope of the Products



Synthetic results



Contact

Contact

Serhiy V. Ryabukhin, Prof. Dr. Sci.
s.v.ryabukhin@gmail.com

Dmytro M. Volochnyuk, Prof. Dr. Sci.
d.volochnyuk@gmail.com

Enamine Ltd, www.enamine.net

78 Winston Churchill St., 02094 Kyiv, Ukraine

References

- ## References
